




Single Coated Foam Tape



- Backing Carrier:** 1/16" white, crosslinked polyethylene foam
- Adhesive:** Highly Moisture-Resistant, Extended-wear Acrylic (M110)
- Release Liner:** 74lb bleached polycoated kraft, silicone treated on one side

Applications

- Electrodes: ECG, EKG, TENS
- Medical device mounting
- AED defibrillator pads
- Transdermal delivery system

Features

Benefits

Polyethylene Foam Backing

- Soft and flexible for maximum patient comfort
- Conformable and easily adapts to medical devices and body contours

Acrylic Adhesive

- Aggressive adhesive designed for extended wear applications
- Removes cleanly without causing trauma
- Highly resistant to perspiration-borne moisture

74lb Kraft Liner

- Easily die-cut into various shapes and sizes
- Good dimensional stability and lay flat

Typical Test Values

	Standard UOM	Metric UOM	Test Method
Total Thickness (with Liner)	70.1 mils	1.781 mm	PSTC 133
Carrier Thickness	62.5 mils	1.588 mm	PSTC 133
Adhesive Thickness	2 mils	0.051 mm	PSTC 133
Liner Thickness	5.6 mils	0.142 mm	PSTC 133
Adhesion to Steel	80 oz/in	8.75 N/cm	PSTC 101 Mod
Recommended Storage Temperature	60 to 80 °F	15.6 to 26.7 °C	
Recommended Storage Humidity	40 to 60 %		