




Double-Coated Closed-Cell Foam Tape



- Backing Carrier:** 1/16" white, closed cell crosslinked polyethylene foam
- Adhesive AS/LS:** Non-sensitizing acrylic (M109)
- Release Liner:** 60lb bleached super calendered kraft, silicone release treated on both sides

Applications

- Electrodes: ECG, EKG, TENS
- Medical device assembly
- Medical device attachment

Features

Polyethylene Foam Carrier

Non-Sensitizing Acrylic Adhesive

60lb Release Liner

Benefits

- Soft and flexible for maximum patient comfort
- Conformable to adapt to medical devices and body contours
- High tack and good shear for a variety of applications
- Stable to EtO and gamma sterilization
- Easy to die-cut and convert

Typical Test Values

	Standard UOM	Metric UOM	Test Method
Total Thickness (with Liner)	68.5 mils	1.74 mm	PSTC 133
Carrier Thickness	62.5 mils	1.588 mm	PSTC 133
Adhesive Thickness (Air Side)	1.1 mils	0.028 mm	PSTC 133
Adhesive Thickness (Liner Side)	1.4 mils	0.036 mm	PSTC 133
Liner Thickness	3.5 mils	0.089 mm	PSTC 133
Adhesion to Steel (Air Side)	160 oz/in	17.5 N/cm	PSTC-101
Adhesion to Steel (Liner Side)	160 oz/in	17.5 N/cm	PSTC 101 Mod
Recommended Storage Temperature	60 to 80 °F	15.6 to 26.7 °C	
Recommended Storage Humidity	40 to 60 %		

Values shown are based on laboratory testing and represent typical product characteristics. Actual performance will vary depending on application conditions, and no guarantee of performance is made. The user is responsible for determining suitability for the intended application. No express or implied warranty, including merchantability or fitness for a particular purpose, is made.